

Stylus™

A COMPLETE LINE OF

OPHTHALMIC KNIVES



Stylus™ Knives

- Trusted In Ophthalmology
- Proprietary Finishing Technology
- Precision Manufacturing

SURGICAL
SPECIALTIES
CORPORATION™

Stab Knives

Stylus™ stab knives are heat-tempered to provide a strong, true point for easy penetration during initial stab incisions.



Knife

52-1501
52-2201
52-3001
52-4501
52-1502
52-2202

Description

15°, Full handle
22.5°, Full handle
30°, Full handle
45°, Full handle
15°, Blade and tip
22.5°, Blade and tip



52-1531
52-1551

15°, Full handle, 3.0mm restricted blade depth
15°, Full handle, 5.0mm restricted blade depth



52-1532
52-1552

15°, Blade and tip, 3.0mm restricted blade depth
15°, Blade and tip, 5.0mm restricted blade depth

Precision Depth Knives

Stylus™ precision depth knives are designed to create groove and limbal relaxation incisions. Single formed design ensures precision accuracy.



Knife


52-3503

Description

350µm, Straight, full handle

Clear Corneal Knives

Stylus™ Clear Corneal knives have been designed specifically for clear corneal incisions. Depth markings of 1.5 mm, 1.75 mm, and 2.0 mm provide easy visual reference for near-clear and hinge incisions.



Knife


52-2663
52-2763
52-2863
52-3061

Description

2.65mm, Angled, full handle with depth markings
2.75mm, Angled, full handle with depth markings
2.85mm, Angled, full handle with depth markings
3.0mm, Angled, full handle with depth markings

Clear Corneal Knives, Double Bevel Parallel Sides

Stylus™ Clear Corneal knives with parallel sides have been designed specifically for clear corneal incisions. Depth markings of 1.5 mm, 1.75 mm, and 2.0 mm provide easy visual reference for near-clear and hinge incisions.



Knife

52-2735
52-2835
52-3035
52-3235

Description

2.75mm, Angled, full handle with depth markings
2.85mm, Angled, full handle with depth markings
3.0mm, Angled, full handle with depth markings
3.2mm, Angled, full handle with depth markings

Cleartrap™ Trapezoid - Clear Corneal/Implant Knives

Stylus™ ClearTrap™ trapezoid clear corneal/implant knives have been designed specifically for clear corneal incisions. Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.



Knife

55-1214
55-2022
55-2724

Description

1.2 - 1.4mm, Angled, full handle with depth markings
2.0 - 2.2mm, Angled, full handle with depth markings
2.75 - 3.2mm, Angled, full handle with depth markings

Stylus™

Slit Knives

Stylus™ Precision Slit Knives Penetrate Easily While Making An Accurate Width Incision For Optimal Fit To The Phacoemulsification Tip.



Knife

52-2801
52-2901

Description

2.85mm, Straight, full handle
3.0mm, Straight, full handle

Slit Knives, Bevel-Up

Stylus™ angled, bevel up knives penetrate easily while creating an accurate width incision for optimal fit to the phacoemulsification tip.



Knife

52-1861
52-2261
52-2461
52-2561
52-2661
52-2761
52-2861
52-2961
52-3261

Description

1.8mm, Angled, full handle, bevel up
2.2mm, Angled, full handle, bevel up
2.4mm, Angled, full handle, bevel up
2.5mm, Angled, full handle, bevel up
2.65mm, Angled, full handle, bevel up
2.75mm, Angled, full handle, bevel up
2.85mm, Angled, full handle, bevel up
3.0mm, Angled, full handle, bevel up
3.2mm, Angled, full handle, bevel up

52-2662

2.65mm, Angled, blade and tip, bevel up

Slit Knives, Double Bevel

Stylus™ double bevel slit knives allow straight, smooth penetration without tissue drag.



Knife

52-1831
52-2231
52-2431
52-2531
52-2631
52-2831

Description

1.8mm, Angled, full handle, double bevel
2.2mm, Angled, full handle, double bevel
2.4mm, Angled, full handle, double bevel
2.5mm, Angled, full handle, double bevel
2.65mm, Angled, full handle, double bevel
2.8mm, Angled, full handle, double bevel

Slit Knives, Parallel Sides

Stylus™ slit knives with parallel sides are designed to hold tissue at the desired width while entering and exiting the wound.



Knife

52-3212

Description

3.2mm, Angled, full handle, bevel up

Sharptome™ Crescent Design Knives

Stylus™ Sharptome™ crescent knives provide consistent sharpness, which nearly eliminates tissue drag. Ideal knives for creating a smooth scleral shelf to promote a self-sealing wound.



Knife

54-5000

Description

2.0mm, Straight, blade and tip

54-1010

2.0mm, Angled, full handle, bevel up

54-1031

2.0mm, Angled, full handle, double bevel

54-6910

1.25mm, Mini-Crescent, angled, full handle

Stylus™

Spoon Knives

Stylus™ unique 320° cutting surface and ultra-thin blade enhance precision cutting capabilities where multi-directional techniques are used.



Knife

54-2261

Description

2.2mm, Angled, full handle, double bevel, 320° cutting edge

54-3061

3.0mm, Angled, full handle, double bevel, 320° cutting edge

Vitrectomy Knives

Stylus™ vitrectomy knives create small, precise access incisions which help prevent leakage around probes and other posterior segment instruments.



Knife

51-2001

Description

20g, Full handle



51-2061

20g, Angled, full handle

Specialty Knives

Stylus™ specialty knives are available in a variety of shapes and sizes to meet the unique needs of various techniques and procedures.



Knife

52-1001

Description

1.0mm, Paracentesis knife, straight, full handle



58-6900

Micro-blade, straight, blade and tip

Guarded Knives

Stylus™ guarded knife features our top-quality Stylus™ blade with the security of a sliding protective sheath. They are available in a wide variety of blade styles including paracentesis, slit, clear corneal and crescent designs.



Knife	Description
52-1501G	15° Stab, Straight, guarded
52-2763G	2.75mm, Clear Corneal, angled, double-bevel, guarded
52-2863G	2.85mm, Clear Corneal, angled, double-bevel, guarded
52-2231G	2.2mm, Slit angled, double-bevel, guarded
52-2431G	2.4mm, Slit angled, double-bevel, guarded
52-2561G	2.5mm, Slit angled, bevel-up, guarded
52-2661G	2.65mm, Slit angled, bevel-up, guarded
52-2761G	2.75mm, Slit angled, bevel-up, guarded
52-2861G	2.85mm, Slit angled, bevel-up, guarded
52-2961G	3.0mm, Slit angled, bevel-up, guarded
54-1010G	2.0mm, Crescent, angled, bevel-up, guarded
58-2010G	1.1mm ClearPort™ paracentesis, angled, guarded

Stylus™



Visit SurgicalSpecialties.com For More Information

Surgical Specialties Corporation, 247 Station Drive, Suite NE1, Westwood, MA 02090, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are owned by the company and its affiliates. ©2001-2016 Surgical Specialties Corporation.
All Rights Reserved. SSC-043 R3 02/18

SURGICAL
SPECIALTIES
CORPORATION™